

THE PROCESS OF A FULL APPLICATION CYCLE:

A whole course of applications* is necessary to exploit the full potential of the EQUUSIR BEST-BOX. This consists of three appointments, although valuable findings can be obtained even after the first one. Each appointment includes comprehensive advice and a detailed report.

Appointment 1:

Initial application, duration: approx. 30 min | Price: 115 €

Each course of applications (application cycle) starts with the initial application. This takes two minutes.

The objective is to take a snapshot of current energy levels: this identifies current interference fields, highlights hidden problem areas, and shows all parts of the body that are working efficiently and are well supplied with energy.



Appointment 2:

Control application, duration: approx. 30 min | Price: 89 €

This determines whether your dog's health has changed since the first application, or whether the measures taken had the intended effect.



Appointment 3:

Final application, duration: approx. 30 min | Price: 89 €

The three reports are compared with each other during the third and final application. From changes over time, conclusions can be drawn regarding long-term measures for stabilising your dog's health, as well as any necessary preventative actions and treatment.

For more information and a appointment booking:
Stephanie Schubert | Phone: +49 151-746 22 941
www.equusir.com

* The EQUUSIR BEST-BOX helps the dogs's own energy fields to recover and become harmonized. It is not therapy within the meaning of medical treatment administered by a professional and should not be regarded as a substitute for a vet. The application complements a medical diagnosis and therapy but is not recognized by medical science.

WHAT DO EXPERTS SAY ABOUT EQUUSIR BEST-BOX?

The EQUUSIR BEST-BOX analyses whether the energy flow is disrupted and eliminates interference fields.

Founder, EQUUSIR Research & Development
Gerold Reinwald



We are really impressed by the concept of the EQUUSIR BEST-BOX and can use it in parallel with our dog behaviour therapy.

Dog trainer
Steve Kaye



In the scientific study we carried out in 2017 on the impact of the EQUUSIR BEST-BOX, we were able to show positive influences over the vegetative nervous system.

Vet / Institute of Experimental and Clinical Pharmacology, Medical University of Graz
Mag. med. vet. Georg Racic



EQUUSIR
THE BEST FOR YOUR DOG

HEALTH BY INNOVATION

EQUUSIR BEST-BOX DOG



EQUUSIR BODY ENERGY
SOLUTION &
BEST-BOX TECHNOLOGY
www.equusir.com


EQUUSIR Germany GmbH
www.equusir.com | info@equusir.com | +49 9605 9199 585
Paul-Engel-Straße 1 | 92729 Weiherhammer | Germany

WHY TREATMENT USING THE EQUUSIR BEST-BOX?



SUPPORT OF THE MEDICAL DIAGNOSIS

External symptoms may have invisible causes and this can be indicated by disharmony in the energy fields.



MONITORING APPLICATION MEASURES

The physical effectiveness of measures can be visualized.



UNDERSTANDING BEHAVIOURAL DISORDERS

Nervousness, fear, unrest in certain situations



RELIEVE ACUTE SYMPTOMS

Problems with skin and fur, muscular and joint diseases, the lymphatic system, respiratory problems



DETECT CHRONIC ILLNESS

Allergies, general fluctuations in immunity and vitality, problems with the skin and fur texture, rheumatism



PROACTIVELY IDENTIFY WEAKNESSES

Cycle irregularities, aggression, fear, dominating behaviour, traumas



RESOLVE INTERFERENCE FIELDS

Muscular tension, neurological abnormalities, strengthening of the muscular and motor system, fascias and connective tissues


HOW IT WORKS:

THE BODY-ENERGY-SOLUTION & TECHNOLOGY (BEST)-BOX USES ENERGY INFORMATION SYSTEMS (EIS®) TO WORK WITH THE BODY'S OWN FREQUENCIES.

DETECT – UNDERSTAND – OPTIMIZE.

These are the terms that describe the EQUUSIR BEST-BOX. The aim is to detect possible weaknesses in the body, subsequently understand these, and optimize them with the help of specialists. Two systems are used in the EQUUSIR BEST-BOX: Infrared-B and coloured light, also known as frequency therapy or photobiomodulation (PBM).

The EQUUSIR BEST-BOX combines modern, high-end technology with the latest findings in quantum physics together with Health & Performance Management (HPM).




The objective of the application is to meet the requirements of modern performance management: Maintain physical health in the long term, intensify the relationship between a dog and its owner and increase physical performance.

With the Body-Energy-Solution & Technology, the EQUUSIR BEST-BOX is used to determine the current status of the energy flow within the body. A detailed report shows which parts of the body are well and adequately supplied with energy, and which ones are blocked.

How can possible weaknesses in the organism be detected with the EQUUSIR BEST-BOX? Imagine healthy cells that exchange information between each other. This exchange runs evenly, at a constant frequency. If an interference field occurs in the organism, this exchange of information is weakened. Possible interference fields can be inflammations, muscular tension, back and joint problems, bacteria, a lack of nutrients, as well as mental load such as stress. The originally constant frequency is automatically off balance. The EQUUSIR BEST-BOX can localize precisely these interference fields and show them in a report.

WHAT DOES A TYPICAL APPLICATION REPORT LOOK LIKE?

The values that are determined are displayed both visually and numerically. For example:



Furthermore, specific problem areas including recommended treatments are assigned to the different centres of the body. Here is a sample extract from the reports:

1. root energy center

- strengthen kidneys and urinary tract
- balance and strengthen the skeletal system
- release tension in the left toe extensor area

2. sacral energy center

- optimize energy flow
- resolve temperamental issues
- balance out power play

3. solar plexus energy center

- transform apathetic behavior
- transforming eating as a purpose in life
- dissolve the tendency of worm infestations
- solve problems in the metabolism of the liver

[...]

HEALTH & PERFORMANCE MANAGEMENT (HPM) - HOLISTICALLY WITH ALL EQUUSIR PRODUCTS

All EQUUSIR systems are designed to support the healing process by adjusting the balance of energy. The EQUUSIR BEST-BOX enables interference fields to be identified precisely and resolved gradually; The BIOS system controls the activity level via neurostimulation and the BIONIC care acts as a disinfectant and regulates the acid-base balance.

Specialists regard Health & Performance Management (HPM) from EQUUSIR as the "interface" between health, prevention and performance. Due to its innovative concepts, EQUUSIR can inform the relevant "team member", which includes the practitioner, physio/osteopath, chiropractor, orthopaedist, nutrition expert, dog trainer and many more, to take specific action.