

WHAT

Instantly heals (micro-)muscle cracks and makes horses loose and supple when being tensed and stressed.

- External muscle and tendon recovering therapy for horses
- Recovers and relaxes fatigued muscles and tendons
- Relieves from muscle cramps
- Approved and certified product by FEI laboratory
- Contains **NO prohibited substances**
- Leaves no sticky residue and color shine behind




3 REASONS WHY

- ***Provides a faster recovery during intense training or competition periods***
Our gel stimulates a fast(er) regeneration process of fatigued muscles and tendons, allowing your horse(s) to recover and return to it's normal condition point far more faster than ever before!
- ***Focus on healing minor muscle and tendon damage***
Arnica and horse Chestnut are known as "counter-irritants" ingredients that ao. triggers the blood circulation (and therefore stimulates the recovery process), relieves tensed and fatigued muscles and promotes the repair of light to mid-heavy impacts.
- ***Our gels are full of natural ingrediënts to accelerate the recovery process of animals and humans***
Minuscule muscle cracks are caused by repetitive (over)loads on specific (muscle) area's in ie. legs, neck and back. Which negatively affects the performance your horse. The recovery gel will enhance and accelerate the healing process through direct stimulation of the affected muscle and tendons.

APPROVED & CERTIFIED

- The Equivi gels have gone through extensive doping tests by an official European FEI laboratory and are proven to be **free of prohibited substances**.
- The principal ingrediënts of the Equivi recovering gel are Waterbase (739mg/g), Rosemary (20mg/g), Arnica (10mg/g) and horse Chestnut (9mg/g).
- We do strongly encourage to check/inform with the respective jurors and their veterinarians at competition to the use of any medication.
- We recommend to not use more than 3 times in a 24h period.
(working efficiency of 2 to 3 hours)



PRODUCT FEATURES

- **DOSAGE & DIRECTIONS FOR USE:**
Equivi recovering gel should be used at full strength and applied liberally (a hand full) on the fatigued or tensed area's and rub in gently via massage until a shiny surface appearance ; allow then to withdraw slowly.
- **STORAGE:** Store below or at room temperature in a cool place. Keep the lid tightly closed.
- **EXPIRY:** 12 months after opening.